

## Product datasheet for **TR519832**

### **lqcg Mouse shRNA Plasmid (Locus ID 69707)**

#### **Product data:**

|                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | shRNA Plasmids                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Product Name:             | lqcg Mouse shRNA Plasmid (Locus ID 69707)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Locus ID:                 | 69707                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Synonyms:                 | 2400003L07Rik; esgd12d; G1-374-12; repro1; stubby12d                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Vector:                   | pRS (TR20003)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| E. coli Selection:        | Ampicillin                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Mammalian Cell Selection: | Puromycin                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Format:                   | Retroviral plasmids                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Components:               | lqcg - Mouse, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 69707). 5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| RefSeq:                   | <a href="#">BC049990</a> , <a href="#">NM_178378</a> , <a href="#">NM_178378.1</a> , <a href="#">NM_178378.2</a> , <a href="#">NM_178378.3</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| UniProt ID:               | <a href="#">Q80W32</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Summary:                  | Component of the nexin-dynein regulatory complex (N-DRC), a key regulator of ciliary/flagellar motility which maintains the alignment and integrity of the distal axoneme and regulates microtubule sliding in motile axonemes. Binds calmodulin when cellular Ca(2+) levels are low and thereby contributes to the regulation of calcium and calmodulin-dependent protein kinase IV (CAMK4) activity; contributes to the regulation of CAMK4 signaling cascades (By similarity). Required for normal axoneme assembly in sperm flagella, normal sperm tail formation and for male fertility (PubMed:24362311, PubMed:24849454). [UniProtKB/Swiss-Prot Function] |
| shRNA Design:             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .                                                                                                                                                                                                                                                                                                               |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).