

## Product datasheet for **TR702011**

### Cc2d1a Rat shRNA Plasmid (Locus ID 288908)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Cc2d1a Rat shRNA Plasmid (Locus ID 288908)
Locus ID:	288908
Synonyms:	RGD1306108
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Cc2d1a - Rat, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 288908). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001013869</a> , <a href="#">NM_001013869.1</a> , <a href="#">BC081948</a>
UniProt ID:	<a href="#">Q66HA5</a>
Summary:	Transcription factor that binds specifically to the DRE (dual repressor element) and represses 5-HT1A gene transcription through this element. Mediates HDAC-independent repression of HTR1A promoter. CAMK2G inhibits CC2D1A-induced repression of the HTR1A. May play a role in the altered regulation of 5-HT1A receptors associated with anxiety and major depression. Performs essential function in controlling functional maturation of synapses (By similarity). [UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).