

Product datasheet for **TR709133**

Nlrp1a Rat shRNA Plasmid (Locus ID 360557)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Nlrp1a Rat shRNA Plasmid (Locus ID 360557)
Locus ID:	360557
Synonyms:	Nalrp1; Nlrp1
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Nlrp1a - Rat, 4 unique 29mer shRNA constructs in retroviral untagged vector (Gene ID = 360557). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	NM_001145755 , NM_001145755.1 , NM_001145755.2
UniProt ID:	D9I2F9
Summary:	Acts as the sensor component of the Nlrp1a inflammasome, which mediates inflammasome activation in response to various pathogen-associated signals, leading to subsequent pyroptosis (By similarity). Inflammasomes are supramolecular complexes that assemble in the cytosol in response to pathogens and other damage-associated signals and play critical roles in innate immunity and inflammation (By similarity). Acts as a recognition receptor (PRR): recognizes specific pathogens and other damage-associated signals, such as B.anthraxis lethal toxin (LT) or Val-boroPro inhibitor, and mediates the formation of the inflammasome polymeric complex (PubMed:20502689, PubMed:31383852). In response to pathogen-associated signals, the N-terminal part of Nlrp1a is degraded by the proteasome, releasing the cleaved C-terminal part of the protein (NACHT, LRR and PYD domains-containing protein 1a, C-terminus), which polymerizes to initiate the formation of the inflammasome complex: the inflammasome directly recruits pro-caspase-1 (proCASP1) independently of PYCARD/ASC and promotes caspase-1 (CASP1) activation, which subsequently cleaves and activates inflammatory cytokines IL1B and IL18 and gasdermin-D (GSDMD), leading to pyroptosis (By similarity).[UniProtKB/Swiss-Prot Function]


[View online »](#)

shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .
Performance Guaranteed:	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>