

## Product datasheet for **TR710974**

### Cebpb Rat shRNA Plasmid (Locus ID 24253)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Cebpb Rat shRNA Plasmid (Locus ID 24253)
Locus ID:	24253
Synonyms:	Il6dbp; NF-IL6; TCF5
Vector:	pRS (TR20003)
E. coli Selection:	Ampicillin
Mammalian Cell Selection:	Puromycin
Format:	Retroviral plasmids
Components:	Cebpb - Rat, 4 unique 29mer shRNA constructs in retroviral untagged vector(Gene ID = 24253). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pRS Vector, TR30012, included for free.
RefSeq:	<a href="#">NM_001301715</a> , <a href="#">NM_001301720</a> , <a href="#">NM_024125</a> , <a href="#">NM_024125.1</a> , <a href="#">NM_024125.2</a> , <a href="#">NM_024125.3</a> , <a href="#">NM_024125.4</a> , <a href="#">NM_024125.5</a> , <a href="#">NM_001301720.1</a> , <a href="#">NM_001301715.1</a> , <a href="#">BC129071</a>
UniProt ID:	<a href="#">P21272</a>
Summary:	This intronless gene encodes a member of the transcription factor family whose members contain a basic leucine-zipper domain. The encoded protein functions as a homodimer but can also form heterodimers with CCAAT/enhancer-binding proteins alpha, delta, and gamma. The encoded protein plays important roles in several cellular processes and in various diseases, including regulating cell proliferation, differentiation, apoptosis and neuroinflammation, and being involved in brain injury and inflammatory progression. The use of alternative in-frame AUG start codons results in multiple protein isoforms, each with different cellular localizations and distinct biological functions. [provided by RefSeq, Sep 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).