

## Product datasheet for **TS400035**

### EXOSC3 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of EXOSC3, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC200035
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	EXOSC3 Mouse Monoclonal Antibody [Clone ID: OTI1D6] (TA504794)
ACCN:	<a href="#">NM_016042</a> , <a href="#">NP_057126</a>
Synonyms:	bA3J10.7; CGI-102; hRrp-40; p10; PCH1B; RRP40; Rrp40p
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_057126</a>
Locus ID:	51010
Cytogenetics:	9p13.2
Protein Families:	Stem cell - Pluripotency
Protein Pathways:	RNA degradation



[View online »](#)