

Product datasheet for **TS400120P5**

CNDP2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CNDP2, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC200120
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CNDP2 Mouse Monoclonal Antibody [Clone ID: OT16D3] (TA503631)
ACCN:	NM_018235 , NP_060705
Synonyms:	CN2; CPGL; HEL-S-13; HsT2298; PEPA
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_060705
Locus ID:	55748
Cytogenetics:	18q22.3
Protein Families:	Protease



[View online »](#)