

Product datasheet for **TS400368P5**

eIF2 alpha (EIF2S1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of EIF2S1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC200368
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	eIF2 alpha (EIF2S1) Mouse Monoclonal Antibody [Clone ID: OTI3H4] (TA501313)
ACCN:	NM_004094 , NP_004085
Synonyms:	EIF-2; EIF-2A; EIF-2alpha; EIF2; EIF2A
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_004085
Locus ID:	1965
Cytogenetics:	14q23.3



[View online »](#)