

## Product datasheet for **TS400386P5**

### Neuraminidase (NEU1) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NEU1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC200386
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	Neuraminidase (NEU1) Mouse Monoclonal Antibody [Clone ID: OTI3D4] (TA801727)
ACCN:	<a href="#">NM_000434</a> , <a href="#">NP_000425</a>
Synonyms:	NANH; NEU; SIAL1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_000425</a>
Locus ID:	4758
Cytogenetics:	6p21.33
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Lysosome, Other glycan degradation, Sphingolipid metabolism



[View online »](#)