

Product datasheet for **TS400415P5**

Sorbitol Dehydrogenase (SORD) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SORD, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC200415
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Sorbitol Dehydrogenase (SORD) Mouse Monoclonal Antibody [Clone ID: OTI9E7] (TA500701)
ACCN:	NM_003104 , NP_003095
Synonyms:	HEL-S-95n; RDH; SDH; SORD1; SORDD; XDH
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_003095
Locus ID:	6652
Cytogenetics:	15q21.1
Protein Families:	Druggable Genome
Protein Pathways:	Fructose and mannose metabolism, Metabolic pathways



[View online »](#)