

## Product datasheet for TS400426P5

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **RFC4 CytoSection**

**Product data:** 

**Product Type:** CytoSections

Description: Transient overexpression of RFC4, transcript variant 1, in HEK293T cells, FFPE control for IHC,

ICC and ISH staining, 25 slides per pack

Species: Human **Expression Host:** HEK293T

**Expression cDNA Clone** 

or AA Sequence:

TrueORF Clone RC200426

C-MYC/DDK Tag:

**Detection Antibodies:** DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)

**Target Detection** 

Antibodies:

RFC4 Mouse Monoclonal Antibody [Clone ID: OTI1A8] (TA506183)

ACCN: NM 002916, NP 002907

A1; RFC37 Synonyms:

Storage: Room Temperature

Stability: Slides are guaranteed for a year from the date of receipt if proper storage instructions were

followed.

Preparation: HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells

were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin &

dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut

and mounted on positively charged SuperFrost slides.

This product is for research use only and is not approved for use in humans or in clinical Note:

diagnosis.

RefSeq: NP 002907

Locus ID: 5984

Cytogenetics: 3q27.3

**Protein Families:** Druggable Genome, Stem cell - Pluripotency

**Protein Pathways:** DNA replication, Mismatch repair, Nucleotide excision repair

