

## Product datasheet for **TS400460P5**

### Hydroxysteroid (17 beta) Dehydrogenase 4 (HSD17B4) CytoSection

#### Product data:

<b>Product Type:</b>	CytoSections
<b>Description:</b>	Transient overexpression of HSD17B4, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	TrueORF Clone RC200460
<b>Tag:</b>	C-MYC/DDK
<b>Detection Antibodies:</b>	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
<b>Target Detection Antibodies:</b>	Hydroxysteroid (17 beta) Dehydrogenase 4 (HSD17B4) Mouse Monoclonal Antibody [Clone ID: OTI4C4] (TA507028)
<b>ACCN:</b>	<u><a href="#">NM_000414</a></u> , <u><a href="#">NP_000405</a></u>
<b>Synonyms:</b>	DBP; MFE-2; MFP-2; MPF-2; PRLTS1; SDR8C1
<b>Storage:</b>	Room Temperature
<b>Stability:</b>	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
<b>Preparation:</b>	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
<b>Note:</b>	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
<b>RefSeq:</b>	<u><a href="#">NP_000405</a></u>
<b>Locus ID:</b>	3295
<b>Cytogenetics:</b>	5q23.1
<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	Metabolic pathways, Primary bile acid biosynthesis



[View online >](#)