

Product datasheet for **TS400500P5**

beta 1 Sodium Potassium ATPase (ATP1B1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ATP1B1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC200500
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<u>NM_001677</u> , <u>NP_001668</u>
Synonyms:	ATP1B
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001668</u>
Locus ID:	481
Cytogenetics:	1q24.2
Protein Families:	Transmembrane
Protein Pathways:	Cardiac muscle contraction



[View online »](#)