

## Product datasheet for **TS400548P5**


### SNRNP40 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SNRNP40 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC200548
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<u>NM_004814</u> , <u>NP_004805</u>
Synonyms:	40K; HPRP8BP; PRP8BP; PRPF8BP; SPF38; WDR57
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_004805</u>
Locus ID:	9410
Cytogenetics:	1p35.2
Protein Pathways:	Spliceosome



[View online »](#)

**Product images:**

Immunohistology staining on SNRNP40 overexpressed cytosection TS400548P5 with mouse anti DDK clone OTI11C3 C/N [TA180144] at 1:400 dilution 4C O/N. Antibody staining was achieved with HIER Citrate pH6 , Polink1 with DAB chromogen detection kit ([D11-18]). Positive stain shown with the brown chromogen present. HEK293T cells were transfected with cDNA clone [RC200548], 5 micron sections, 40x magnification