

Product datasheet for **TS400604P5**

Cytochrome b5 Outer Mitochondrial Membrane (CYB5B) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CYB5B in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC200604
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_030579</u> , <u>NP_085056</u>
Synonyms:	CYB5-M; CYPB5M; OMB5
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_085056</u>
Locus ID:	80777
Cytogenetics:	16q22.1
Protein Families:	Transmembrane



[View online »](#)