

Product datasheet for **TS400614**

FH CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of FH in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC200614
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	FH Mouse Monoclonal Antibody [Clone ID: OT11F10] (TA500675)
ACCN:	NM_000143 , NP_000134
Synonyms:	FMRD; HLRCC; HsFH; LRCC; MCL; MCUL1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_000134
Locus ID:	2271
Cytogenetics:	1q43
Protein Families:	Druggable Genome
Protein Pathways:	Citrate cycle (TCA cycle), Metabolic pathways, Pathways in cancer, Renal cell carcinoma



[View online »](#)