

Product datasheet for **TS400663**

Melanoma gp100 (PMEL) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PMEL, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC200663
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Melanoma gp100 (PMEL) Mouse Monoclonal Antibody [Clone ID: OTI10E5] (TA500424)
ACCN:	NM_006928 , NP_008859
Synonyms:	D12S53E; gp100; ME20; ME20-M; ME20M; P1; P100; PMEL17; SI; SIL; SILV
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_008859
Locus ID:	6490
Cytogenetics:	12q13.2
Protein Families:	Secreted Protein, Transmembrane



[View online »](#)