

## Product datasheet for **TS400855P5**

### POLR2F CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of POLR2F in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC200855
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	POLR2F Mouse Monoclonal Antibody [Clone ID: OTI2F2] (TA807758)
ACCN:	<a href="#">NM_021974</a> , <a href="#">NP_068809</a>
Synonyms:	HRBP14.4; POLRF; RPABC2; RPABC14.4; RPB6; RPB14.4; RPC15
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_068809</a>
Locus ID:	5435
Cytogenetics:	22q13.1
Protein Families:	Transcription Factors
Protein Pathways:	Huntington's disease, Metabolic pathways, Purine metabolism, Pyrimidine metabolism, RNA polymerase



[View online »](#)