

Product datasheet for TS400963P5

GRHPR CytoSection

Product data:

Product Type: CytoSections **Description:** Transient overexpression of GRHPR in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack Species: Human **Expression Host:** HEK293T **Expression cDNA Clone** TrueORF Clone RC200963 or AA Sequence: C-MYC/DDK Tag: **Detection Antibodies:** Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100) **Target Detection** GRHPR Mouse Monoclonal Antibody [Clone ID: OTI2H4] (TA502091) Antibodies: ACCN: NM 012203, NP 036335 GLXR; GLYD; PH2 Synonyms: Storage: **Room Temperature** Stability: Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. **Preparation:** HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. This product is for research use only and is not approved for use in humans or in clinical Note: diagnosis. **RefSeq:** NP 036335 Locus ID: 9380 **Cytogenetics:** 9p13.2 **Protein Families:** Druggable Genome **Protein Pathways:** Glyoxylate and dicarboxylate metabolism, Metabolic pathways, Pyruvate metabolism



/iew online >

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn