

Product datasheet for **TS400995P5**

Mutarotase (GALM) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of GALM in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC200995
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	Mutarotase (GALM) Mouse Monoclonal Antibody [Clone ID: OTI2F8] (TA812316)
ACCN:	NM_138801 , NP_620156
Synonyms:	BLOCK25; GALAC4; GLAT; HEL-S-63p; IBD1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_620156
Locus ID:	130589
Cytogenetics:	2p22.1
Protein Pathways:	Glycolysis / Gluconeogenesis



[View online »](#)