

## Product datasheet for **TS401015P5**

### RBMS1 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RBMS1, transcript variant 3, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201015
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	RBMS1 Mouse Monoclonal Antibody [Clone ID: OT12H1] (TA808580)
ACCN:	<a href="#">NM_002897</a> , <a href="#">NP_002888</a>
Synonyms:	C2orf12; HCC-4; MSSP; MSSP-1; MSSP-2; MSSP-3; SCR2; YC1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_002888</a>
Locus ID:	5937
Cytogenetics:	2q24.2



[View online »](#)