

Product datasheet for **TS401025P5**

ZDHHC9 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ZDHHC9, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201025
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_016032 , NP_057116
Synonyms:	CGI89; CXorf11; DHHC9; MMSA1; MRXSZ; ZDHHC10; ZNF379; ZNF380
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_057116
Locus ID:	51114
Cytogenetics:	Xq26.1
Protein Families:	Transmembrane



[View online »](#)