

Product datasheet for **TS401074P5**

TRIM44 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TRIM44 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201074
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	TRIM44 Mouse Monoclonal Antibody [Clone ID: OTI1E10] (TA505561)
ACCN:	NM_017583 , NP_060053
Synonyms:	AN3; DIPB; HSA249128; MC7
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_060053
Locus ID:	54765
Cytogenetics:	11p13



[View online »](#)