

Product datasheet for **TS401225**

MMP7 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MMP7 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201225
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	MMP7 Mouse Monoclonal Antibody [Clone ID: OTI3G7] (TA804842)
ACCN:	NM_002423 , NP_002414
Synonyms:	MMP-7; MPLS1; PUMP-1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_002414
Locus ID:	4316
Cytogenetics:	11q22.2
Protein Families:	Druggable Genome, Protease
Protein Pathways:	Wnt signaling pathway



[View online »](#)