

Product datasheet for **TS401235P5**

ISG15 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ISG15 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201235
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	ISG15 Mouse Monoclonal Antibody [Clone ID: OT16C8] (TA807536)
ACCN:	NM_005101 , NP_005092
Synonyms:	G1P2; hUCRP; IFI15; IMD38; IP17; UCRP
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_005092
Locus ID:	9636
Cytogenetics:	1p36.33
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	RIG-I-like receptor signaling pathway



[View online »](#)