

Product datasheet for **TS401238**

POLD1 CytoSection

Product data:

| | |
|--|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of POLD1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC201238 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| ACCN: | <u>NM_002691</u> , <u>NP_002682</u> |
| Synonyms: | CDC2; CRCS10; MDPL; POLD |
| Storage: | Room Temperature |
| Stability: | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | <u>NP_002682</u> |
| Locus ID: | 5424 |
| Cytogenetics: | 19q13.33 |
| Protein Families: | Druggable Genome, Stem cell - Pluripotency |
| Protein Pathways: | Base excision repair, DNA replication, Homologous recombination, Metabolic pathways, Mismatch repair, Nucleotide excision repair, Purine metabolism, Pyrimidine metabolism |



[View online »](#)