

## Product datasheet for **TS401266**

### POLR2E CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of POLR2E in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201266
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	POLR2E Mouse Monoclonal Antibody [Clone ID: OTI3B5] (TA502548)
ACCN:	<a href="#">NM_002695</a> , <a href="#">NP_002686</a>
Synonyms:	hRPB25; hsRPB5; RPABC1; RPB5; XAP4
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_002686</a>
Locus ID:	5434
Cytogenetics:	19p13.3
Protein Families:	Transcription Factors
Protein Pathways:	Huntington's disease, Metabolic pathways, Purine metabolism, Pyrimidine metabolism, RNA polymerase



[View online »](#)