

Product datasheet for **TS401377**

AP2M1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of AP2M1, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201377
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	AP2M1 Mouse Monoclonal Antibody [Clone ID: OT11E9] (TA503019)
ACCN:	NM_001025205 , NP_001020376
Synonyms:	AP50; CLAPM1; MRD60; mu2
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001020376
Locus ID:	1173
Cytogenetics:	3q27.1
Protein Families:	Druggable Genome
Protein Pathways:	Endocytosis, Huntington's disease



[View online »](#)