

Product datasheet for **TS401397**

Grp75 (HSPA9) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HSPA9 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201397
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Grp75 (HSPA9) Mouse Monoclonal Antibody [Clone ID: OTI9F8] (TA500516)
ACCN:	NM_004134 , NP_004125
Synonyms:	CRP40; CSA; EVPLS; GRP-75; GRP75; HEL-S-124m; HSPA9B; MOT; MOT2; MTHSP75; PBP74; SAAN; SIDBA4
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_004125
Locus ID:	3313
Cytogenetics:	5q31.2
Protein Families:	Stem cell - Pluripotency
Protein Pathways:	RNA degradation



[View online »](#)