

## Product datasheet for **TS401601**

### SULT1A1 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SULT1A1, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201601
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	SULT1A1 Mouse Monoclonal Antibody [Clone ID: OTI5F1] (TA501837)
ACCN:	<a href="#">NM_001055</a> , <a href="#">NP_001046</a>
Synonyms:	HAST1/HAST2; P-PST; P-PST 1; PST; ST1A1; ST1A3; STP; STP1; ts-PST; TSPST1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001046</a>
Locus ID:	6817
Cytogenetics:	16p11.2
Protein Pathways:	Sulfur metabolism



[View online »](#)