

## Product datasheet for **TS401663P5**

### CEP72 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CEP72 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201663
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	CEP72 Mouse Monoclonal Antibody [Clone ID: OTI4D2] (TA811145)
ACCN:	<u><a href="#">NM_018140</a></u> , <u><a href="#">NP_060610</a></u>
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_060610</a></u>
Locus ID:	55722
Cytogenetics:	5p15.33



[View online »](#)