

Product datasheet for **TS401725P5**

HYPE (FICD) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of FICD in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201725
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	HYPE (FICD) Mouse Monoclonal Antibody [Clone ID: OTI2A11] (TA802841)
ACCN:	NM_007076 , NP_009007
Synonyms:	HIP13; HYPE; UNQ3041
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_009007
Locus ID:	11153
Cytogenetics:	12q23.3
Protein Families:	Transmembrane



[View online »](#)