

Product datasheet for **TS401738P5**

Thimet Oligopeptidase (THOP1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of THOP1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201738
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	Thimet Oligopeptidase (THOP1) Mouse Monoclonal Antibody [Clone ID: OTI2D3] (TA501409)
ACCN:	NM_003249 , NP_003240
Synonyms:	EP24.15; MEPD_HUMAN; MP78; TOP
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_003240
Locus ID:	7064
Cytogenetics:	19p13.3
Protein Families:	Druggable Genome, Protease
Protein Pathways:	Renin-angiotensin system



[View online »](#)