

Product datasheet for **TS401803**

PGP9.5 (UCHL1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of UCHL1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201803
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	PGP9.5 (UCHL1) Mouse Monoclonal Antibody [Clone ID: OTI2B5] (TA504289)
ACCN:	NM_004181 , NP_004172
Synonyms:	HEL-117; HEL-S-53; NDGOA; PARK5; PGP 9.5; PGP9.5; PGP95; SPG79; Uch-L1
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_004172
Locus ID:	7345
Cytogenetics:	4p13
Protein Families:	Druggable Genome, Protease
Protein Pathways:	Parkinson's disease



[View online »](#)