

## Product datasheet for **TS401813**

### PEX16 Cytosection

#### Product data:

Product Type:	Cytosections
Description:	Transient overexpression of PEX16 (NM_057174) in HEK293T cells paraffin embedded 4 um sections, controls for ICC/IHC staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201813
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Accession Number:	<u>NM_057174</u> , <u>NP_476515</u>
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 4 um sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_476515</u>
Locus ID:	9409
Cytogenetics:	11p11.2



[View online »](#)