

Product datasheet for **TS401847P5**

UAP56 (DDX39B) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DDX39B, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201847
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	UAP56 (DDX39B) Mouse Monoclonal Antibody [Clone ID: OTI2C5] (TA501186)
ACCN:	NM_080598 , NP_542165
Synonyms:	BAT1; D6S81E; UAP56
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_542165
Locus ID:	7919
Cytogenetics:	6p21.33
Protein Pathways:	Spliceosome



[View online »](#)