

Product datasheet for **TS401945P5**

PPP2R5D CytoSection

Product data:

| | |
|---------------------------------------|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of PPP2R5D, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC201945 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| Target Detection Antibodies: | PPP2R5D Mouse Monoclonal Antibody [Clone ID: OTI5E7] (TA803916) |
| ACCN: | NM_006245 , NP_006236 |
| Synonyms: | B56D; B56delta; MRD35 |
| Storage: | Room Temperature |
| Stability: | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | NP_006236 |
| Locus ID: | 5528 |
| Cytogenetics: | 6p21.1 |
| Protein Families: | Phosphatase |
| Protein Pathways: | Oocyte meiosis, Wnt signaling pathway |



[View online »](#)