

Product datasheet for **TS401974**

Carbonic Anhydrase II (CA2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CA2 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC201974
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_000067 , NP_000058
Synonyms:	CA-II; CAC; CAII; Car2; HEL-76; HEL-S-282
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_000058
Locus ID:	760
Cytogenetics:	8q21.2
Protein Families:	Druggable Genome
Protein Pathways:	Nitrogen metabolism



[View online »](#)