

## Product datasheet for **TS402064**

### SIAT4A (ST3GAL1) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ST3GAL1, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202064
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u><a href="#">NM_173344</a></u> , <u><a href="#">NP_775479</a></u>
Synonyms:	1; Gal-NAc6S; SIAT4A; SIATFL; ST3GalA; ST3GalA.1; ST3GalIA; ST3O
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_775479</a></u>
Locus ID:	6482
Cytogenetics:	8q24.22
Protein Families:	Secreted Protein, Transmembrane
Protein Pathways:	Glycosphingolipid biosynthesis - ganglio series, Glycosphingolipid biosynthesis - globo series, Keratan sulfate biosynthesis, Metabolic pathways, O-Glycan biosynthesis



[View online »](#)