

Product datasheet for **TS402167**

GCIP interacting protein p29 (SYF2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SYF2, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202167
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	GCIP interacting protein p29 (SYF2) Mouse Monoclonal Antibody [Clone ID: OT16A5] (TA808453)
ACCN:	NM_015484 , NP_056299
Synonyms:	CBPIN; fSAP29; NTC31; P29
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_056299
Locus ID:	25949
Cytogenetics:	1p36.11
Protein Pathways:	Spliceosome



[View online »](#)