

Product datasheet for **TS402172**

Cytochrome P450 Reductase (POR) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of POR in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202172
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Cytochrome P450 Reductase (POR) Mouse Monoclonal Antibody [Clone ID: OT16C4] (TA500601)
ACCN:	NM_000941 , NP_000932
Synonyms:	CPR; CYPOR; P450R
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_000932
Locus ID:	5447
Cytogenetics:	7q11.23
Protein Families:	Druggable Genome, P450, Transmembrane



[View online »](#)