

Product datasheet for **TS402222**

MAP1LC3A CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MAP1LC3A, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202222
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	MAP1LC3A Mouse Monoclonal Antibody [Clone ID: OT110H6] (TA803577)
ACCN:	NM_032514 , NP_115903
Synonyms:	ATG8E; LC3; LC3A; MAP1ALC3; MAP1BLC3
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_115903
Locus ID:	84557
Cytogenetics:	20q11.22



[View online »](#)