

Product datasheet for **TS402227**

LCMT2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of LCMT2 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202227
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_014793 , NP_055608
Synonyms:	PPM2; TYW4
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_055608
Locus ID:	9836
Cytogenetics:	15q15.3
Protein Families:	Druggable Genome
Protein Pathways:	Alzheimer's disease, Androgen and estrogen metabolism, Cardiac muscle contraction, Histidine metabolism, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease, Selenoamino acid metabolism, Tyrosine metabolism



[View online »](#)