

## Product datasheet for **TS402309**

### GAPDH CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of GAPDH in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202309
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	GAPDH Mouse Monoclonal Antibody [Clone ID: OTI2D9] (TA802519)
ACCN:	<u><a href="#">NM_002046</a></u> , <u><a href="#">NP_002037</a></u>
Synonyms:	G3PD; GAPD; HEL-S-162eP
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_002037</a></u>
Locus ID:	2597
Cytogenetics:	12p13.31
Protein Families:	ES Cell Differentiation/IPS
Protein Pathways:	Alzheimer's disease, Glycolysis / Gluconeogenesis, Metabolic pathways



[View online »](#)