

## Product datasheet for **TS402333P5**

### ATF 4 (ATF4) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ATF4, transcript variant 2, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202333
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	ATF 4 (ATF4) Mouse Monoclonal Antibody [Clone ID: OT11H9] (TA812797)
ACCN:	<a href="#">NM_182810</a> , <a href="#">NP_877962</a>
Synonyms:	CREB-2; CREB2; TAXREB67; TXREB
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_877962</a>
Locus ID:	468
Cytogenetics:	22q13.1
Protein Families:	Transcription Factors
Protein Pathways:	GnRH signaling pathway, Long-term potentiation, MAPK signaling pathway, Neurotrophin signaling pathway, Prostate cancer



[View online »](#)