

## Product datasheet for **TS402348**

### CANT1 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CANT1, transcript variant 1, in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 5 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202348
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CANT1 Mouse Monoclonal Antibody [Clone ID: OT17F8] (TA805629)
ACCN:	<a href="#">NM_138793</a> , <a href="#">NP_620148</a>
Synonyms:	DBQD; DBQD1; EDM7; SCAN-1; SCAN1; SHAPY
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_620148</a>
Locus ID:	124583
Cytogenetics:	17q25.3
Protein Families:	Secreted Protein
Protein Pathways:	Purine metabolism, Pyrimidine metabolism



[View online »](#)