

Product datasheet for **TS402365P5**

SLUG (SNAI2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SNAI2 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202365
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	SLUG (SNAI2) Mouse Monoclonal Antibody [Clone ID: OTI1A6] (TA800167)
ACCN:	NM_003068 , NP_003059
Synonyms:	SLUG; SLUGH; SLUGH1; SNAI2; WS2D
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_003059
Locus ID:	6591
Cytogenetics:	8q11.21
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Adherens junction



[View online »](#)