

Product datasheet for **TS402837P5**

acyl CoA Thioesterase 2 (ACOT2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ACOT2 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202837
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_006821 , NP_006812
Synonyms:	CTE-IA; CTE1A; MTE1; PTE2; PTE2A; ZAP128
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 μ m sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_006812
Locus ID:	10965
Cytogenetics:	14q24.3
Protein Pathways:	Biosynthesis of unsaturated fatty acids



[View online »](#)