

## Product datasheet for **TS402885**

### LKB1 (STK11) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of STK11 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC202885
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	LKB1 (STK11) Mouse Monoclonal Antibody [Clone ID: OTI10A3] (TA802554)
ACCN:	<a href="#">NM_000455</a> , <a href="#">NP_000446</a>
Synonyms:	hLKB1; LKB1; PJS
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_000446</a>
Locus ID:	6794
Cytogenetics:	19p13.3
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Adipocytokine signaling pathway, mTOR signaling pathway



[View online »](#)