

Product datasheet for **TS403097P5**

H4C5 CytoSection

Product data:

| | |
|---------------------------------------|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of H4C5 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC203097 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| ACCN: | <u>NM_003545</u> , <u>NP_003536</u> |
| Synonyms: | H4-16; H4/j; H4C1; H4C2; H4C3; H4C4; H4C6; H4C8; H4C9; H4C11; H4C12; H4C13; H4C14; H4C15; H4Fj; HIST1H4E |
| Storage: | Room Temperature |
| Stability: | Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | <u>NP_003536</u> |
| Locus ID: | 8367 |
| Cytogenetics: | 6p22.2 |
| Protein Pathways: | Systemic lupus erythematosus |



[View online »](#)