

Product datasheet for **TS403261P5**

C15orf15 (RSL24D1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RSL24D1 in HEK293T cells, FFPE control for IHC, ICC and ISH staining, 25 slides per pack
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC203261
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<u>NM_016304</u> , <u>NP_057388</u>
Synonyms:	C15orf15; HRP-L30-iso; L30; RLP24; RPL24; RPL24L; TVAS3
Storage:	Room Temperature
Stability:	Slides are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. 5 µm sections of the FFPE cell pellet blocks are cut and mounted on positively charged SuperFrost slides.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_057388</u>
Locus ID:	51187
Cytogenetics:	15q21.3
Protein Pathways:	Ribosome



[View online »](#)